


ORDERING GUIDE

	CODE	DESCRIPTION	NO./PACK
SINGLE-USE	EU3783	SPRINK S.H.F. 125 mm GRID, complete KIT	10
	EU3784	SPRINK S.H.F. 160 mm GRID, complete KIT	10
	EU3785	SPRINK S.H.F. 195 mm GRID, complete KIT	10
	EU5234	Sprink LED Light	10
REUSABLE	EU3235	Pressure reducer for Sprink	1
	Wall socket for connection to the compressed air source		
	EU1026	UNI9507 compressed air connector	1
	EU1156	SIO connector (AFNOR NR-S-90-116)	1
	EU1157	DIN1326 compressed air connector	1
	EU1158	BS5682 compressed air connector	1



SPRINK SINGLE-USE KIT

- 1 Lavage ON/OFF clamp
- 2 Upper lavage cannula
- 3 Lower lavage cannula
- 4 Splash shield
- 5 Suction hose
- 6 Saline bag spike connection
- 7 Pressure regulator Luer-lock connection
- 8 Sterilizing air filter
- 9 Sprink LED Light

