

CANNULAE



CANNULAE PORTFOLIO

FOR ADULT
& PAEDIATRIC
PATIENTS



EUROSETS

EVERY LIFE MATTERS

TABLE OF CONTENTS

ARTERIAL CANNULAE	3
Arterial Cannulae Wire-Reinforced	4
Aortic Catheters	6
VENOUS CANNULAE	7
Single Stage Venous Cannulae	8
Dual Stage Venous Cannulae	9
Triple Stage Venous Cannulae	9
CARDIOPLEGIA CANNULAE	10
Aortic Root Cannulae	11
Ostial Perfusion Cannulae	13
Retrograde Cardioplegia Cannulae	15
Cardioplegia Accessories	16
VENTS & SUCTION DEVICES	19
Vent Catheters	20
Suction Devices	24
MICS CANNULAE	27
Arterial Cannulae	28
Venous Cannulae	31
Cardioplegia Cannulae	32
ECLS CANNULAE	33
Venous and Dual Lumen Cannulae	34
Arterial Cannulae	38
Femoral Perfusion Set	39
PAEDIATRIC CANNULAE	40
Arterial Cannulae	41
Venous Cannulae	44
ACCESSORIES	46
Tourniquets	47
Connectors	48
Vessel Cannulae	50
Gas Diffusers	50
Insertion Kits	51



CANNULAE



ARTERIAL CANNULAE


- ♥ ARTERIAL CANNULAE WIRE-REINFORCED
- ♥ AORTIC CATHETERS



Arterial Cannulae Wire-Reinforced


Eurosets Arterial Cannulae are intended for returning oxygenated blood to the patient. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for arterial direct cannulation.

Arterial Cannulae, Wire-Reinforced, Straight tip with sewing collar



SIZE (FR)	WITHOUT CONNECTOR	WITH 3/8" LUER CONNECTOR AND VENT CAP	LENGTH	COATING	PCS/BOX
14	FLKA22A14	FLKA22D14	230 mm	NO	10
16	FLKA22A16	FLKA22D16	230 mm	NO	10
18	FLKA22A18	FLKA22D18	230 mm	NO	10
20	FLKA22A20	FLKA22D20	230 mm	NO	10
22	FLKA22A22	FLKA22D22	230 mm	NO	10
24	FLKA22A24	FLKA22D24	230 mm	NO	10
26	FLKA22A26	FLKA22D26	230 mm	NO	10


Arterial Cannulae, Wire-Reinforced, Straight tip with ring



SIZE (FR)	WITH 3/8" LUER CONNECTOR AND VENT CAP	LENGTH	COATING	PCS/BOX
14	FLKA24D14	230 mm	NO	10
16	FLKA24D16	230 mm	NO	10
18	FLKA24D18	230 mm	NO	10
20	FLKA24D20	230 mm	NO	10
22	FLKA24D22	230 mm	NO	10
24	FLKA24D24	230 mm	NO	10
26	FLKA24D26	230 mm	NO	10




Arterial Cannulae, Wire-Reinforced, Curved tip with ring



SIZE (FR)	WITH 3/8" LUER CONNECTOR AND VENT CAP	LENGTH	COATING	PCS/BOX
14	FLKA09D14	230 mm	NO	10
16	FLKA09D16	230 mm	NO	10
18	FLKA09D18	230 mm	NO	10
20	FLKA09D20	230 mm	NO	10
22	FLKA09D22	230 mm	NO	10
24	FLKA09D24	230 mm	NO	10
26	FLKA09D26	230 mm	NO	10

Arterial Cannulae, Wire-Reinforced, Curved tip with sewing collar




SIZE (FR)	WITHOUT CONNECTOR	WITH 3/8" LUER CONNECTOR AND VENT CAP	LENGTH	COATING	PCS/BOX
14	FLKA07A14	FLKA07D14	230 mm	NO	10
16	FLKA07A16	FLKA07D16	230 mm	NO	10
18	FLKA07A18	FLKA07D18	230 mm	NO	10
20	FLKA07A20	FLKA07D20	230 mm	NO	10
22	FLKA07A22	FLKA07D22	230 mm	NO	10
24	FLKA07A24	FLKA07D24	230 mm	NO	10
26	FLKA07A26	FLKA07D26	230 mm	NO	10



Aortic Catheters

Eurosets Non-Reinforced Aortic Catheters are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass. Available in a variety of configurations, these cannulae features an easy tip orientation thanks to the indexing marker and a smooth bevelled tip.

Aortic Perfusion Cannula, Non-Reinforced, Straight hard tip with ring



SIZE (FR)	CODE	CONNECTION	LENGTH	COATING	PCS/BOX
24	FLKA32D24	3/8" LUER LOCK	200 mm	NO	10
24	FLKA36D24	3/8" LUER LOCK	300 mm	NO	10
24	FLKA34D24	3/8" LUER LOCK	400 mm	NO	10



CANNULAE



VENOUS CANNULAE


- ▼ SINGLE STAGE VENOUS CANNULAE
- ▼ DUAL STAGE VENOUS CANNULAE
- ▼ TRIPLE STAGE VENOUS CANNULAE



Single Stage Venous Cannulae


Eurosets Venous Cannulae are intended for the drainage of deoxygenated blood from the patient. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for venous direct cannulation.

Single Stage Venous Cannulae, Wire-Reinforced, Bullet Tip



SIZE (FR)	STRAIGHT TIP REINFORCED	STRAIGHT TIP REINFORCED WITH RING	90° CURVED TIP REINFORCED	90° CURVED TIP REINFORCED WITH FLEXIBLE RING	90° CURVED SHORT TIP WITH RING	LENGTH	CONNECTION SITE	COATING	PCS/BOX
20	FLKV02A20	FLKV04A20	FLKV01A20	FLKV05A20	FLKV03A20	370 mm	1/4"	NO	10
22	FLKV02A22	FLKV04A22	FLKV01A22	FLKV05A22	FLKV03A22	370 mm	1/4"	NO	10
24	FLKV02A24	FLKV04A24	FLKV01A24	FLKV05A24	FLKV03A24	370 mm	1/4"	NO	10
26	FLKV02A26	FLKV04A26	FLKV01A26	FLKV05A26	FLKV03A26	370 mm	3/8"	NO	10
28	FLKV02A28	FLKV04A28	FLKV01A28	FLKV05A28	FLKV03A28	370 mm	3/8"	NO	10
30	FLKV02A30	FLKV04A30	FLKV01A30	FLKV05A30	FLKV03A30	370 mm	3/8"	NO	10
32	FLKV02A32	FLKV04A32	FLKV01A32	FLKV05A32	FLKV03A32	370 mm	3/8"	NO	10
34	FLKV02A34	FLKV04A34	FLKV01A34	FLKV05A34	FLKV03A34	370 mm	3/8"	NO	10
36	FLKV02A36	FLKV04A36	FLKV01A36	FLKV05A36	FLKV03A36	370 mm	3/8"	NO	10
38	FLKV02A38	FLKV04A38	FLKV01A38	FLKV05A38	FLKV03A38	370 mm	3/8"	NO	10
40	FLKV02A40	FLKV04A40	FLKV01A40	FLKV05A40	FLKV03A40	370 mm	3/8"	NO	10

Single Stage Venous Cannulae, Wire-Reinforced, Curved metal tip





SIZE (FR)	CURVED METAL TIP	LENGTH	CONNECTION SITE	COATING	PCS/BOX
20	FLKV07A20	388 mm	1/4"	NO	10
22	FLKV07A22	388 mm	1/4"	NO	10
24	FLKV07A24	398 mm	1/4"	NO	10
26	FLKV07A26	410 mm	3/8"	NO	10
28	FLKV07A28	410 mm	3/8"	NO	10
30	FLKV07A30	420 mm	3/8"	NO	10
32	FLKV07A32	420 mm	3/8"	NO	10



Dual Stage Venous Cannulae

Eurosets Venous Dual Stage Cannulae are intended for the drainage of deoxygenated blood from the Right Atrium & IVC. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for venous direct cannulation.

Dual Stage Venous Cannulae





SIZE (FR)	ROUND BODY WITHOUT CONNECTOR	ROUND BODY WITH CONNECTOR	FLAT BODY WITHOUT CONNECTOR	FLAT BODY WITH CONNECTOR	CONNECTION SITE	LENGTH	COATING	PCS/BOX
28/36	FLKV06A2836	FLKV06F2836	FLKV12A2836	FLKV12F2836	1/2"	370 mm	NO	10
30/38	FLKV06A3038	FLKV06F3038	FLKV12A3038	FLKV12F3038	1/2"	370 mm	NO	10
32/40	FLKV06A3240	FLKV06F3240	FLKV12A3240	FLKV12F3240	1/2"	380 mm	NO	10
34/46	FLKV06A3446	FLKV06F3446	FLKV12A3446	FLKV12F3446	1/2"	380 mm	NO	10
36/46	FLKV06A3646	FLKV06F3646	FLKV12A3646	FLKV12F3646	1/2"	390 mm	NO	10
36/51	FLKV06A3651	FLKV06F3651	FLKV12A3651	FLKV12F3651	1/2"	390 mm	NO	10
38/46	FLKV06A3846	FLKV06F3846	FLKV12A3846	FLKV12F3846	1/2"	390 mm	NO	10

Triple Stage Venous Cannulae

Eurosets Venous Triple Stage Cannulae are intended for the drainage of deoxygenated blood during VAVD procedures. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for venous direct cannulation.

Triple Stage Venous Cannulae



SIZE (FR)	ROUND BODY WITHOUT CONNECTOR	ROUND BODY WITH CONNECTOR	FLAT BODY WITHOUT CONNECTOR	FLAT BODY WITH CONNECTOR	CONNECTION SITE	LENGTH	COATING	PCS/BOX
28/36/36	FLKV08A2836	FLKV08F2836	FLKV14A2836	FLKV14F2836	1/2"	380 mm	NO	10
32/40/40	FLKV08A3240	FLKV08F3240	FLKV14A3240	FLKV14F3240	1/2"	400 mm	NO	10
34/46/46	FLKV08A3446	FLKV08F3446	FLKV14A3446	FLKV14F3446	1/2"	400 mm	NO	10
36/46/46	FLKV08A3646	FLKV08F3646	FLKV14A3646	FLKV14F3646	1/2"	400 mm	NO	10



CANNULAE



CARDIOPLEGIA CANNULAE


- ▼ AORTIC ROOT CANNULAE
- ▼ OSTIAL PERFUSION CANNULAE
- ▼ RETROGRADE CARDIOPLEGIA CANNULAE
- ▼ CARDIOPLEGIA ACCESSORIES



Aortic Root Cannulae


Eurosets Aortic Root Cannulae are used during cardiopulmonary bypass for antegrade delivery of cardioplegia. An optional vent line allows to perform venting of the heart when needed. A variety of configurations and sizes are available to suit different patients and types of operations.

Aortic Root Cannulae



NEEDLE SIZE (GA)	WITHOUT VENT LINE	WITH VENT LINE	TIP SIZE (FR)	LENGTH	PCS/BOX
16	EUARS70007	EUARS80007	6	150 mm	20
14	EUARS70008	EUARS80008	7	150 mm	20
12	EUARS70010	EUARS80010	9	150 mm	20

Aortic Root Cannulae with anti-bleeding Valve



NEEDLE SIZE (GA)	WITHOUT VENT LINE	WITH VENT LINE	TIP SIZE (FR)	LENGTH	PCS/BOX
16	EUARS70007V	EUARS80007V	6	150 mm	20
14	EUARS70008V	EUARS80008V	7	150 mm	20
12	EUARS70010V	EUARS80010V	9	150 mm	20



Aortic Root Cannulae with Blue Tip



NEEDLE SIZE (GA)	WITHOUT VENT LINE	WITH VENT LINE AND WHITE CLAMP	WITH VENT LINE AND RED AND BLUE CLAMPS	WITH TOP VENT LINE AND WHITE CLAMP	WITH TOP VENT LINE AND RED AND BLUE CLAMPS	LENGTH	PCS/BOX
14	FLKC20K16	FLKC04K16	-	FLKC12K16	FLKC16K16	170 mm	20
12	FLKC20K20	FLKC04K20	FLKC08K20	FLKC12K20	FLKC16K20	170 mm	20
11	-	FLKC04K23	FLKC08K23	-	-	170 mm	20

Aortic Root Cannulae "alternative version"




NEEDLE SIZE (GA)	WITHOUT VENT LINE	LENGTH WITHOUT VENT LINE	WITH VENT LINE	LENGTH WITH VENT LINE	PCS/BOX
18	FLKC54K10	170 mm	FLKC52K10	155 mm	20
16	FLKC54K13	170 mm	FLKC52K13	155 mm	20
14	FLKC54K16	170 mm	FLKC52K16	155 mm	20
12	FLKC54K20	170 mm	FLKC52K20	155 mm	20



Ostial Perfusion Cannulae


Eurosets Ostial Perfusion Cannulae are intended to delivery cardioplegia directly into the coronary arteries during cardiopulmonary bypass. Available with different tips, sizes and with self-inflating balloon versions, Eurosets Ostial Perfusion Cannulae can suit every cardioplegia delivery need.

Ostial Perfusion Cannulae with Basket Tip



SIZE (FR)	45° TIP	90° TIP	LENGTH	CONNECTION	PCS/BOX
10	EUOP10010	EUOP20010	210 mm	FEMALE LUER LOCK	20
12	EUOP10012	EUOP20012	210 mm	FEMALE LUER LOCK	20
14	EUOP10014	EUOP20014	210 mm	FEMALE LUER LOCK	20

Ostial Perfusion Cannulae with Self-Inflating Balloon Tip and Soft Body



BALLOON SIZE (mm)	TIP SIZE	CURVED	LENGTH	CONNECTION	PCS/BOX
6	3 mm	EUOP9006SF	300 mm	FEMALE LUER LOCK	20
7	3 mm	EUOP9007SF	300 mm	FEMALE LUER LOCK	20



Ostial Perfusion Cannulae with Soft Tip





SIZE (mm)	135° TIP	90° TIP	LENGTH	CONNECTION	PCS/BOX
3	FLKC03130	FLKC01130	237 mm	FEMALE LUER LOCK	10
3,5	FLKC03135	FLKC01135	237 mm	FEMALE LUER LOCK	10
4	FLKC03140	FLKC01140	237 mm	FEMALE LUER LOCK	10



Retrograde Cardioplegia Cannulae

Eurosets Retrograde Cardioplegia Cannulae are intended to be connected to the cardioplegia line in order to infuse retrogradely cardioplegic solution and blood into the patient's coronary sinus during CPB procedure.

Retrograde Cardioplegia Cannulae with Self Inflating Balloon and Guidewire Stylet



SIZE (FR)	CODE	BALLOON DIAMETER	LENGTH	PCS/BOX
14	EURC40014G	18 mm	310 mm	20



Cardioplegia Accessories

Eurosets has a wide range of Cardioplegia Accessories and adapters intended to aid in the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass surgery.













Cardioplegia Needle












NEEDLE	CODE	SIZE (GA)	LENGTH	CONNECTION	PCS/BOX
14 mm	EUCN00016	16	50 mm	LUER LOCK	20



Cardioplegia Adapters

CODE	DESCRIPTION		PCS/BOX
EUCAL00003	Perfusion adapter with 2 male luer lock and 1 female luer lock and color-coded clamps		20
EUCAL00004	Multiple perfusion adapter for aortic root and three vein grafts		20
EUCAL00005	Multiple perfusion adapter for aortic root and three vein grafts with attached vessel cannulae		20
EUCAL00006	Dual perfusion set with female luer lock, cone connector and color coded clamps		20
EUCAL00007	Multiple perfusion adapter with vent line and vessel cannula		20
EUCAL00008	Aortic perfusion set with vent (long = 20 cm)		20
EUCAL00009	Dual perfusion set with female luer lock, male luer lock connector, color coded clamps		20
EUCAL00010	Multiple set for perfusion of three vein grafts		20
EUCAL00011	Multiple set for perfusion of three vein grafts with attached vessel cannulae		20
EUCAL00012	Multiple perfusion adapter for retrograde and aortic root with vent		20
EUCAL00013	Aortic perfusion set with vent (short = 11 cm)		20
EUCAL00014	Aortic perfusion set with 1/4x1/4x1/4 connector		20



CODE	DESCRIPTION	PCS/BOX
EUCAL00015	Straight adapter, 14 cm 	20
EUCAL00016	Straight adapter, 15 cm 	20
EUCAL00017	Perfusion adapter 1/4" 	20
EUCAL00018	Perfusion adapter 3/8" 	20
EUCAL00019	Perfusion adapter 1/4" with 3-way stopcock 	20
EUCAL00020	Perfusion adapter with 2 female luer lock and 1 male luer lock and color-coded clamps 	20
EUCAL00021	Recirculation y-adapter 	20
EUCAL00022	Multiple perfusion adapter with vent line 	20
EUCAL00023	Six multiple perfusion adapter 	20



CANNULAE



VENTS & SUCTION DEVICES

- ▼ VENT CATHETERS
- ▼ SUCTION DEVICES



Vent Catheters

Eurosets Vent Catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.

Silicone Vent Catheter with Closed Tip



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
16	EUSV40116G	65 mm	350 mm	20
20	EUSV40120G	65 mm	350 mm	20

PVC Vent Catheters, Straight, Open Tip, Without Stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
17	EUVCS90317	60 mm	320 mm	20

PVC Vent Catheters, Straight, Open Tip, With plastic stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
18	EUVCS90318R	60 mm	320 mm	20



PVC Vent Catheters, Straight, Open Tip, With guidewire stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
18	EUVCS90318G	60 mm	320 mm	20

PVC Vent Catheters, Straight, Closed Tip, Without stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
18	EUVCS40318	60 mm	320 mm	20

PVC Vent Catheters, Straight, Closed, With plastic stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
18	EUVCS40318R	60 mm	320 mm	20
20	EUVCS40320R	60 mm	320 mm	20



PVC Vent Catheters, Straight, Closed Tip, With guidewire stylet



SIZE (FR)	CODE	TIP PERFORATION	LENGTH	PCS/BOX
18	EUVCS40318G	60 mm	320 mm	20
20	EUVCS40320G	60 mm	320 mm	20

Left Vent Catheter with guidewire



SIZE (FR)	CODE	CONNECTION SITE	LENGTH	PCS/BOX
12	FLKX00212	1/4"	350 mm	20
14	FLKX00214	1/4"	350 mm	20
17	FLKX00217	1/4"	350 mm	20

Left Vent Catheter with PVC introducer



SIZE (FR)	CODE	CONNECTION SITE	LENGTH	PCS/BOX
12	FLKX00412	1/4"	350 mm	20
17	FLKX00417	1/4"	350 mm	20



Left Vent Catheter, With stainless steel tube on the tip



SIZE (FR)	CODE	LENGTH	CONNECTION SITE	PCS/BOX
12	FLKX01412	270 mm	1/4"	20
17	FLKX01417	320 mm	1/4"	20

PVC Vent Catheters, Angled, Without stylet



SIZE (FR)	CODE	TIP PERFORATION	CONNECTION SITE	LENGTH	PCS/BOX
13	EUVCA40313	30 mm	1/4"	320 mm	20
20	EUVCA40320	100 mm	1/4"	340 mm	20

Left Vent Cannula, Angled tip



SIZE (FR)	CODE	CONNECTION SITE	LENGTH	PCS/BOX
10	FLKX01310	1/4"	310 mm	20
13	FLKX01313	1/4"	310 mm	20



Suction Devices

Eurosets Suction Devices are medical devices designed to remove fluid from the surgical field during cardiopulmonary bypass procedures.

Suction Wands



CODE	TIP	LENGTH	CONNECTOR	PCS/BOX
EURS100300	FLUTED TIP	360 mm	1/4"	20

BOIDA, Mini & Rigid Suckers



CODE	TYPE	LENGTH	CONNECTOR	PCS/BOX
FLKS01C20	BOIDA	290 mm	1/4"	20
FLKS02C20	RIGID	290 mm	1/4"	20

Intracardiac suckers



SIZE (FR)	CODE	TUBING	LENGTH	CONNECTOR	PCS/BOX
17	EUFS100317	FLEXIBLE	360 mm	1/4"	20



Pericardial Sumps with metal tip



SIZE (FR)	CODE	TUBING	TIP LENGTH	LENGTH	CONNECTOR	PCS/BOX
18	EUPS120316	FLEXIBLE	20 mm	300 mm	1/4"	20
23	EUPS130316	FLEXIBLE	35 mm	300 mm	1/4"	20

Pericardial Sumps with wirewound tip



SIZE (FR)	CODE	TUBING	TIP LENGTH	LENGTH	CONNECTOR	PCS/BOX
16	EUPS110316	FLEXIBLE	40 mm	300 mm	1/4"	20



Pericardial Sump with protective spring



SIZE (FR)	CODE	LENGTH	CONNECTOR	PCS/BOX
17	FLKX01617	675 mm	1/4"	20



CANNULAE



MICS CANNULAE



- ▼ ARTERIAL CANNULAE
- ▼ VENOUS CANNULAE
- ▼ SMART DRAINAGE DEVICE
- ▼ CARDIOPLEGIA CANNULAE



Arterial Cannulae

Eurosets MICS Arterial Cannulae are intended for returning oxygenated blood to the patient through the femoral artery or the ascending aorta. Optimized for a percutaneous or direct insertion, they are available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, representing the optimal choice for arterial cannulation in minimally invasive procedures.



Femoral Arterial Cannulae with flexible ring



SIZE (FR)	CODE	CONNECTOR	LENGTH	COATING	PCS/BOX
12	FLKA10D12	3/8" LUER LOCK	235 mm	NO	5
14	FLKA10D14	3/8" LUER LOCK	235 mm	NO	5
16	FLKA10D16	3/8" LUER LOCK	235 mm	NO	5
18	FLKA10D18	3/8" LUER LOCK	235 mm	NO	5
20	FLKA10D20	3/8" LUER LOCK	235 mm	NO	5
22	FLKA10D22	3/8" LUER LOCK	235 mm	NO	5
24	FLKA10D24	3/8" LUER LOCK	235 mm	NO	5



Femoral Arterial Cannulae “Long Version”



SIZE (FR)	CODE	CONNECTOR	LENGTH	COATING	PCS/BOX
14	FLKA08D14	3/8" LUER LOCK	285 mm	NO	5
16	FLKA08D16	3/8" LUER LOCK	285 mm	NO	5
18	FLKA08D18	3/8" LUER LOCK	285 mm	NO	5
20	FLKA08D20	3/8" LUER LOCK	285 mm	NO	5
22	FLKA08D22	3/8" LUER LOCK	285 mm	NO	5
24	FLKA08D24	3/8" LUER LOCK	285 mm	NO	5



“Elonga” Arterial Cannulae with straight tip and flexible ring





SIZE (FR)	CODE	CONNECTOR	LENGTH	COATING	PCS/BOX
14	FLKA06D14	3/8" LUER LOCK	235 mm	NO	5
16	FLKA06D16	3/8" LUER LOCK	235 mm	NO	5
18	FLKA06D18	3/8" LUER LOCK	235 mm	NO	5
20	FLKA06D20	3/8" LUER LOCK	235 mm	NO	5
22	FLKA06D22	3/8" LUER LOCK	235 mm	NO	5
24	FLKA06D24	3/8" LUER LOCK	235 mm	NO	5



Venous Cannulae

Eurosets Venous Cannulae are intended for femoral cannulation to drain blood from the right atrium and/or vena cava during cardiopulmonary bypass procedures. They are characterized by an excellent profile in terms of biocompatibility, superior performances and easy handling for a smoother insertion and positioning.

Venous Femoral Cannulae




SIZE (FR)	CODE	CONNECTOR	LENGTH	COATING	PCS/BOX
20	FLKV16B20	3/8"	750 mm	NO	1
22	FLKV16B22	3/8"	750 mm	NO	1
24	FLKV16B24	3/8"	750 mm	NO	1
26	FLKV16B26	3/8"	750 mm	NO	1
28	FLKV16B28	3/8"	750 mm	NO	1



Cardioplegia Cannulae

A dedicated offering of long aortic root cannulae, without vent line, is available to deliver cardioplegia from remote accesses.

Aortic Root Long Needle



NEEDLE SIZE (GA)	CODE	TIP SIZE (FR)	VENT LINE	LENGTH	COATING	PCS/BOX
14	EUARL70008	7	NO	280 mm	NO	20



CANNULAE



ECLS CANNULAE


- ♥ VENOUS AND DUAL LUMEN CANNULAE
- ♥ ARTERIAL CANNULAE
- ♥ FEMORAL PERFUSION SET



Venous and Dual Lumen Cannulae

Eurosets Venous Ecmo Cannulae are intended for the drainage of deoxygenated blood from the patient. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for the venous peripheral cannulation. Validated up to 30 days.



Revas Veno Venous Dual Lumen Cannulae



SIZE (FR)	CODE	CONNECTION SIZE	LENGTH	COATING	PCS/BOX
19	FLKW19R	1/4"	305 mm	YES	1
23	FLKW23R	3/8"	480 mm	YES	1
27	FLKW27R	3/8"	480 mm	YES	1
31	FLKW31R	3/8"	480 mm	YES	1





Revas Universal Ecmo/Ecls Cannulae



SIZE (FR)	CODE	CONNECTION SITE	LENGTH	COATING	PCS/BOX
18	FLKV20B18R	3/8"	930 mm	YES	1
20	FLKV20B20R	3/8"	930 mm	YES	1
22	FLKV20B22R	3/8"	930 mm	YES	1





Revas Venous Femoral Cannulae



CONFIGURATION	CODE	SIZE (FR)	CONNECTION SITE	LENGTH	COATING	PCS/BOX
SHORT VERSION	FLKV10SB22R	22	3/8"	625 mm	YES	1
	FLKV10SB24R	24	3/8"	625 mm	YES	1
	FLKV10SB26R	26	3/8"	625 mm	YES	1
LONG VERSION	FLKV10B22R	22	3/8"	820 mm	YES	1
	FLKV10B24R	24	3/8"	930 mm	YES	1
	FLKV10B26R	26	3/8"	930 mm	YES	1



Revas Ecmo/Ecls Venous Femoral Cannulae Pre Bent-Tip





SIZE (FR)	CODE	CONNECTION SITE	LENGTH	COATING	PCS/BOX
18	FLKV19B18R	3/8"	86 cm	YES	1
20	FLKV19B20R	3/8"	86 cm	YES	1
22	FLKV19B22R	3/8"	86 cm	YES	1



Arterial Cannulae

Eurosets Arterial Ecmo Cannulae are intended to reinfuse oxygenated blood to the patient. Available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, they represent the optimal choice for the arterial peripheral cannulation. Validated up to 30 days.

Revas Arterial Femoral Cannulae





CONFIGURATION	SIZE (FR)	CODE	CONNECTION SIZE	LENGTH	COATING	PCS/BOX
SHORT VERSION	16	FLKA18D16R	3/8" LL	265 mm	YES	5
	18	FLKA18D18R	3/8" LL	265 mm	YES	5
	20	FLKA18D20R	3/8" LL	265 mm	YES	5
	22	FLKA18D22R	3/8" LL	265 mm	YES	5
	24	FLKA18D24R	3/8" LL	265 mm	YES	5
LONG VERSION	16	FLKA18LD16R	3/8" LL	355 mm	YES	5
	18	FLKA18LD18R	3/8" LL	355 mm	YES	5
	20	FLKA18LD20R	3/8" LL	355 mm	YES	5
	22	FLKA18LD22R	3/8" LL	355 mm	YES	5
	24	FLKA18LD24R	3/8" LL	355 mm	YES	5



Femoral Perfusion Set

Eurosets Femoral Perfusion Set is intended for the leg-perfusion during ECMO procedures involving a cannulation of the femoral artery.

Crurasave - Femoral Perfusion Set



SIZE (FR)	CODE	CONFIGURATION	COATING	PCS/BOX
7	FLKX00107	Cannula with introducer, guidewire, 18 GA introducer needle, syringe	NO	5



CANNULAE



PAEDIATRIC CANNULAE


- ▼ ARTERIAL CANNULAE
- ▼ VENOUS CANNULAE



Arterial Cannulae


Eurosets Paediatric Arterial Cannulae are intended for returning oxygenated blood to the patient through the femoral artery or the ascending aorta. Optimized for the smaller patients, they are available in a variety of configurations, kink-resistant, with low pressure drop value and DEHP-Free tubing, representing the optimal choice for arterial cannulation in paediatric CPB procedures.

Arterial Cannulae, Wire-Reinforced



SIZE (FR)	WITH 3/16" CONNECTOR AND LUER LOCK	WITH 1/4" CONNECTOR AND LUER LOCK	LENGTH	COATING	PCS/BOX
8	FLKA12G08	FLKA12E08	200 mm	NO	5
10	FLKA12G10	FLKA12E10	200 mm	NO	5
12	FLKA12G12	FLKA12E12	200 mm	NO	5
14	FLKA12G14	FLKA12E14	200 mm	NO	5
16	FLKA12G16	FLKA12E16	200 mm	NO	5
18	FLKA12G18	FLKA12E18	200 mm	NO	5


Femoral Arterial Cannulae, Wire-Reinforced, With flexible ring



SIZE (FR)	WITH 1/4" CONNECTOR AND LUER LOCK	LENGTH	COATING	PCS/BOX
8	FLKA10E08	140 mm	NO	5
10	FLKA10E10	140 mm	NO	5
12	FLKA10E12	140 mm	NO	5
14	-	140 mm	NO	5




Arterial Cannulae, Wire-Reinforced, Straight tip with sewing collar



SIZE (FR)	WITH 1/4" LL CONNECTOR	WITHOUT CONNECTOR	LENGTH	COATING	PCS/BOX
8	FLKA22E08	FLKA22A08	230 mm	NO	10
10	FLKA22E10	FLKA22A10	230 mm	NO	10
12	FLKA22E12	FLKA22A12	230 mm	NO	10
14	FLKA22E14	FLKA22A14	230 mm	NO	10


Arterial Cannulae, Wire-Reinforced, Straight tip with ring



SIZE (FR)	WITH 1/4" LL CONNECTOR	LENGTH	COATING	PCS/BOX
8	FLKA24E08	230 mm	NO	10
10	FLKA24E10	230 mm	NO	10
12	FLKA24E12	230 mm	NO	10
14	FLKA24E14	230 mm	NO	10




Arterial Cannulae, Wire-Reinforced, Curved tip with ring



SIZE (FR)	WITH 1/4" LL CONNECTOR	LENGTH	COATING	PCS/BOX
8	FLKA09E08	230 mm	NO	10
10	FLKA09E10	230 mm	NO	10
12	FLKA09E12	230 mm	NO	10
14	FLKA09E14	230 mm	NO	10

Arterial Cannulae, Wire-Reinforced, Curved tip with sewing collar




SIZE (FR)	WITHOUT CONNECTOR	WITH 1/4" LL CONNECTOR	LENGTH	COATING	PCS/BOX
8	FLKA07A08	FLKA07E08	230 mm	NO	10
10	FLKA07A10	FLKA07E10	230 mm	NO	10
12	FLKA07A12	FLKA07E12	230 mm	NO	10
14	FLKA07A14	FLKA07E14	230 mm	NO	10



Venous Cannulae


Eurosets Paediatric Venous Cannulae are intended to drain deoxygenated blood from the patient's body. Available in a variety of configurations to suits even the smallest vessels, they are characterized by low pressure drop values, kink-resistant tubings and optimal biocompatibility, representing the optimal choice for paediatric CPB procedures.

Venous Femoral Cannulae, Wire-Reinforced



SIZE (FR)	WITH 3/16" CONNECTOR	WITH 1/4" CONNECTOR	LENGTH	COATING	PCS/BOX
8	FLKV18J08	FLKV18C08	200 mm	NO	1
10	FLKV18J10	FLKV18C10	200 mm	NO	1
12	FLKV18J12	FLKV18C12	200 mm	NO	1
14	FLKV18J14	FLKV18C14	200 mm	NO	1


Venous Cannulae, Wire-Reinforced, Straight Tip



SIZE (FR)	WITH SUTURE RING	WITHOUT SUTURE RING	LENGTH	CONNECTION	COATING	PCS/BOX
12	FLKV04A12	FLKV02A12	370 mm	1/4"	NO	10
14	FLKV04A14	FLKV02A14	370 mm	1/4"	NO	10
16	FLKV04A16	FLKV02A16	370 mm	1/4"	NO	10
18	FLKV04A18	FLKV02A18	370 mm	1/4"	NO	10
20	FLKV04A20	FLKV02A20	370 mm	1/4"	NO	10
22	FLKV04A22	FLKV02A22	370 mm	1/4"	NO	10
24	FLKV04A24	FLKV02A24	370 mm	1/4"	NO	10




Venous Cannulae, Wire-Reinforced, Curved Tip



SIZE (FR)	WITHOUT SUTURE RING	WITH FLEXIBLE RING	SHORT TIP WITH RING	LENGTH	CONNECTION	COATING
12	FLKV01A12	FLKV05A12	FLKV03A12	370 mm	1/4"	NO
14	FLKV01A14	FLKV05A14	FLKV03A14	370 mm	1/4"	NO
16	FLKV01A16	FLKV05A16	FLKV03A16	370 mm	1/4"	NO
18	FLKV01A18	FLKV05A18	FLKV03A18	370 mm	1/4"	NO
20	FLKV01A20	FLKV05A20	FLKV03A20	370 mm	1/4"	NO
22	FLKV01A22	FLKV05A22	FLKV03A22	370 mm	1/4"	NO
24	FLKV01A24	FLKV05A24	FLKV03A24	370 mm	1/4"	NO

Venous Cannulae, Wire-Reinforced, Curved Metal Tip



SIZE (FR)	CODE	LENGTH	CONNECTION	COATING	PCS/BOX
12	FLKV07A12	372 mm	1/4"	NO	10
14	FLKV07A14	384 mm	1/4"	NO	10
16	FLKV07A16	380 mm	1/4"	NO	10
18	FLKV07A18	386 mm	1/4"	NO	10
20	FLKV07A20	388 mm	1/4"	NO	10
22	FLKV07A22	388 mm	1/4"	NO	10
24	FLKV07A24	398 mm	1/4"	NO	10



ACCESSORIES

- ▼ TOURNIQUETS
- ▼ CONNECTORS
- ▼ VESSEL CANNULAE
- ▼ GAS DIFFUSERS
- ▼ INSERTION KITS




Tourniquets

Eurosets Cannulation Tourniquets sets are intended to secure arterial and venous cannulae on the proper anatomical site. Available in two different configurations, they represent an optimal choice for cannulation during CPB procedures.

Tourniquets



A



B

SIZE (FR)	CODE	NUMBER OF TUBES	LENGTH	PCS/BOX
9	EU3206 (Tourquetin A)	2	TUBE LENGTH 130 mm	50
11	EU3824 (Tourquetin B)	5	TUBE LENGTH 170 mm, SNARE LENGTH 250 mm	50



Connectors

Eurosets Connectors are intended for the connection of the cannulae to the extracorporeal bypass system during CPB procedures. Several different configurations are available.

Connectors Adapters



PORT 1	PORT 2	PORT 3	WITH LUER	WITHOUT LUER	PCS/BOX
3/16"	3/16"	3/16"	EU3404	EU3361	100
3/16"	3/16"	1/4"	EU3373	EU3372	100
1/4"	1/4"	1/4"	EU3366	EU3358	100
1/4"	1/4"	3/8"	EU3369	EU3362	100
1/4"	3/8"	1/4"	EU3429	EU3356	100
1/4"	3/8"	3/8"	EU3374	EU3364	100
1/4"	3/8"	1/2"	EU3520	EU3519	100
1/4"	1/2"	1/2"	EU3377	EU3375	100
3/8"	3/8"	1/4"	EU3357	EU3376	100
3/8"	3/8"	3/8"	EU3367	EU3359	100
3/8"	3/8"	1/2"	EU3370	EU3363	100
3/8"	1/2"	1/2"	EU3371	EU3365	100
1/2"	1/2"	1/2"	EU3368	EU3360	100

Connectors Adapters

PORT 1	LUER	LUER LOCK	PCS/BOX
1/4"	EU3405	EU3398	100
3/8"	EU3406	EU3399	100



Straight Connectors Adapters



PORT 1	PORT 2	WITH LUER	WITHOUT LUER	PCS/BOX
3/16"	3/16"	EU3388	EU3381	100
3/16"	1/4"	EU3394	EU3393	100
1/4"	1/4"	EU3385	EU3378	100
1/4"	3/8"	EU3389	EU3382	100
1/4"	1/2"	EU3392	EU3384	100
3/8"	3/8"	EU3386	EU3379	100
3/8"	1/2"	EU3390	EU3383	100
3/8"	5/8"	EU3396	EU3391	100
1/2"	1/2"	EU3387	EU3380	100
1/2"	5/8"	EU3397	EU3395	100



Vessel Cannulae

Eurosets Vessel Cannulae are intended for vein graft and/or for the anterograde cardioplegia administration during CPB procedures.

Transparent Vessel Cannulae



LENGTH	LUER	LUER LOCK	PCS/BOX
5 cm	EUVC90000FL	-	40
6,5 cm	-	EUVB90000FL	40

White Vessel Cannulae



LENGTH	LUER	LUER LOCK	PCS/BOX
5 cm	EUVC90000R	-	40

Gas Diffusers

Eurosets Gas Diffuser is intended for intraoperative wound ventilation with Co₂ inflation to prevent air embolism and bacteriological contamination during CPB procedures.

Malleable Gas Diffuser




SHAFT LENGTH	TUBING LENGTH	CONNECTION SITE	CONFIGURATION	CODE	PCS/BOX
40 cm - MALLEABLE	225 cm	1/4"	standard	EUGS130000	20



Insertion Kits

Eurosets Insertion Kits are intended to perform the peripheral insertion of arterial & venous cannulae during CPB or ECLS procedures.



Arterial

CODE	DESCRIPTION	COMPONENTS	GUIDEWIRE	DILATORS SIZES (FR)	PCS/BOX
EURIK80A	Complete configuration with soft & transparent case	18 Ga introducer needle, scalpel, syringe, guidewire, dilator set	80 cm x 0.035"	8-12-16-18-20-22-24	1

Venous

CODE	DESCRIPTION	COMPONENTS	GUIDEWIRE	DILATORS SIZES (FR)	PCS/BOX
EURIK180V	Complete configuration with soft & transparent case	18 Ga introducer needle, scalpel, syringe, guidewire, dilator set	180 cm x 0.035"	8-12-16-18-20-22-24	1

Venous Large

CODE	DESCRIPTION	COMPONENTS	GUIDEWIRE	DILATORS SIZES (FR)	PCS/BOX
EURIK180VL	Complete configuration with soft & transparent case	18 Ga introducer needle, scalpel, syringe, guidewire, dilator set	180 cm x 0.035"	8-12-16-20-24-26-28-30	1

EUROSETS EUROPE

ITALY

Eurosets S.r.l.
Strada Statale 12, 143
41036 Medolla (MO)
Italy
Ph: +39 0535 660311
Fax: +39 0535 51248
info@eurosets.com

BENELUX

Eurosets Benelux
Winstar Park,
Rue Provinciale, 62
1301 Bierges (Wavre)
Belgium
Ph: +32 10237993
info.nl@eurosets.com
info.be@eurosets.com

FRANCE

Eurosets France Sarl
Immeuble LE FONTENOY
96, Boulevard Vivier Merle
69423 LYON Cedex 03
Ph: +33 962191542
Fax: +33 970066573
info.fr@eurosets.com

GERMANY

Eurosets GmbH
Industriestr. 29
82194 Gröbenzell
Ph: +49 81426671227
Fax: +4981426671229
info.de@eurosets.com

UK

Eurosets Ltd
12 Bridewell Place,
Third Floor East
London, EC4V 6AP
United Kingdom
Ph: 07859154270
info@eurosets.com

EUROSETS ASIA

CHINA

Eurosets (Suzhou) Medical Device Ltd.
Suite 1308-2, The Summit (South Block)
No. 119, Suzhou Avenue West,
Suzhou Industrial Park - 215000, Greater China.
苏州工业园区苏州大道西119
号苏悦广场南楼1308-2, 215000,
更大的中国
info.cn@eurosets.com

EUROSETS SOUTH AMERICA

LATAM

Eurosets Latam Sas
NIT: 901.838.279
Dirección del domicilio principal:
Cl 93 No. 11
A - 28 Of 601
Municipio: Bogotá D.C.
info.latam@eurosets.com

EUROSETS EMIRATES

UAE

Eurosets Fzco
2003A, Single Business Tower
Business Bay, Dubai, UAE

   [eurosets.com](https://www.eurosets.com)

Prodotto da
Manufactured by
Fabriqué par
Hersteller
Fabricado por

For "EU" codes:
Andocor n.v.
Kwikaard 104
B-2980 Zoersel, Belgium
Email: info@andocor.com

For "FLK" Codes
Free Life Medical GmbH
Kleinheidstrasse 1a
52080 Aachen - Germany
Tel: +49-2405-419 900-0
Fax: +49-2405-419 900-90
E-mail: info@freelife-gmbh.de

STERILE/EO



CE ₀₁₂₃ Latex free



Distribuito da
Distributed by
Distribué par
Vertrieb
Distribuido por

Eurosets s.r.l.
Strada Statale 12, n°143
41036 Medolla (MO) Italy
Tel: +39 0535 660311
Fax: +39 0535 51248
E-mail: info@eurosets.com
eurosets.com